

## Baricitinib?(Phosphate)

Catalog No: #S2360


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Baricitinib?(Phosphate)
Brief Description	Inhibitors
Purification	98.00%
Target Name	JAK inhibitor; Tyrosine Kinases inhibitor
Calculated MW	469.41
Formulation	C16H20N7O6PS
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:OP(O)(O)=O.CCS(=O)(=O)N1CC(CC#N)(C1)N1C=C(C=N1)C1=C2C=CNC2=NC=N1Pathways:Angiogenesis; Proteases/Proteasome; Chromatin/Epigenetic; JAK/STAT signaling; Stem CellsReceptor:JAK1; JAK2; JAK3; Tyk2Boiling pt:Melting pt:Solubility:200 mM in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Fridman JS, et al. J Immunol, 2010, 184(9), 5298-5307.

Note: This product is for in vitro research use only