

Irbinitinib(ARRY-380; ONT-380)

Catalog No: #S2364


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Irbinitinib(ARRY-380; ONT-380)
Brief Description	Inhibitors
Purification	99.96%
Target Name	EGFR inhibitor; HER inhibitor
Calculated MW	480.52
Formulation	C26H24N8O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CC1=C(OC2=CC3=NC=NN3C=C2)C=CC(NC2=NC=NC3=C2C=C(NC2=NC(C)(C)CO2)C=C3)=C1Pathways:Angiogenesis; JAK/STAT signaling; Tyrosine Kinase/AdaptorsReceptor:EGFR; HER2/ErbB2; p95 HER2Boiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. P Lee, et al. Cancer Research(Poster Session Abstracts). 2009, 69(24 Suppl):Abstract nr 5104.

Note: This product is for in vitro research use only