

Tipiracil hydrochloride

Catalog No: #S2366


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tipiracil hydrochloride
Brief Description	Inhibitors
Purification	99.50%
Target Name	Nucleoside Antimetabolite/Analog inhibitor
Calculated MW	279.12
Formulation	C ₉ H ₁₂ Cl ₂ N ₄ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: Cl.C1=C(CN2CCCC2=N)NC(=O)NC1=O Pathways: DNA Damage/DNA Repair Receptor: Nucleoside Antimetabolite/Analog Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Yoshino T, et al. Lancet Oncol. 2012 Oct;13(10):993-1001.

Note: This product is for in vitro research use only