

## Edoxaban Tosylate Monohydrate

Catalog No: #S2368


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Edoxaban Tosylate Monohydrate
Brief Description	Inhibitors
Purification	99.86%
Target Name	Factor Xa inhibitor; Thrombin inhibitor
Calculated MW	738.27
Formulation	C31H38ClN7O7S2H2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: CC1=CC=C(C=C1)S(=O)(=O)O.CN1CCC2=C(C1)SC(=N2)C(=O)N[C@@H]3C[C@H](CC[C@@H]3NC(=O)C(=O)NC4=NC=C(C=C4)Cl)C(=O)N(C)C

Pathways: Metabolism; Others Receptor: Factor Xa; Factor Xa; Thrombin

Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL (135.45 mM)

Appearance: White Powder

Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Furugohri T, et al. J Thromb Haemost. 2008, 6(9), 1542-1549.

Note: This product is for in vitro research use only