

Emapunil

Catalog No: #S2371


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Emapunil
Brief Description	Inhibitors
Purification	99.80%
Target Name	Others
Calculated MW	401.46
Formulation	C23H23N5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCN(CC1=CC=CC=C1)C(=O)CN1C(=O)N(C)C2=CN=C(N=C12)C1=CC=CC=C1
 Pathways: Others
 Receptor: TSPO
 ligand
 Boling pt: Melting pt: Solubility: DMSO: 21 mg/mL
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Rupprecht R, et al. Science. 2009 Jul 24;325(5939):490-3.

Note: This product is for in vitro research use only