

## GNL3L Polyclonal Antibody

Catalog No: #40969

Package Size: #40969-1 50ul #40969-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	GNL3L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	GNL3L Polyclonal Antibody detects endogenous levels of GNL3L protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GNL3L.
Target Name	GNL3L
Other Names	GNL3L; Guanine nucleotide-binding protein-like 3-like protein
Accession No.	Swiss-Prot: Q9NVN8NCBI Gene ID: 54552
Uniprot	Q9NVN8
GeneID	54552;
SDS-PAGE MW	66kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/5000.

Not yet tested in other applications.

## Images



Western Blot analysis of HT29 cells using GNL3L Polyclonal Antibody



Western Blot analysis of NIH-3T3 LOVO cells using GNL3L Polyclonal Antibody

Note: This product is for in vitro research use only