


## Description

Product Name	LY2835219
Brief Description	Inhibitors
Purification	99.43%
Target Name	CDK inhibitor
Calculated MW	506.59
Formulation	C27H32F2N8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:CCN1CCN(CC2=CN=C(NC3=NC=C(F)C(=N3)C3=CC4=C(N=C(C)N4C(C)C)C(F)=C3)C=C2)CC1Pathways:Cell Cycle/CheckpointReceptor:CDK4; CDK6Boiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Patrick J. Roberts, et al. JNCI J Natl Y Inst (2012) 104 (6): 476-487.

Note: This product is for in vitro research use only