

Ipragliflozin

Catalog No: #S2385


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ipragliflozin
Brief Description	Inhibitors
Purification	99.97%
Target Name	SGLT inhibitor
Calculated MW	404.45
Formulation	C21H21FO5S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: C1=CC=C2C(=C1)C=C(S2)CC3=C(C=CC(=C3)[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O)F Pathways: GPCR/G

Protein Receptor: SGLT2 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Imamura M, et al. Bioorg Med Chem. 2012 May 15;20(10):3263-79.

Note: This product is for in vitro research use only