

GPR101 Polyclonal Antibody

Catalog No: #40970

Package Size: #40970-1 50ul #40970-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR101 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Mk
Specificity	GPR101 Polyclonal Antibody detects endogenous levels of GPR101 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human GPR101.
Target Name	GPR101
Other Names	GPR101; Probable G-protein coupled receptor 101
Accession No.	Swiss-Prot: Q96P66NCBI Gene ID: 83550
Uniprot	Q96P66
GeneID	83550;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using GPR101 Polyclonal Antibody

Note: This product is for in vitro research use only