

Glycopyrrolate

Catalog No: #S2387


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycopyrrolate
Brief Description	Inhibitors
Purification	99.41%
Target Name	AChR antagonist
Calculated MW	398.33
Formulation	C ₁₉ H ₂₈ BrNO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: [Br-].[C][N+](C)CCC(C1)OC(=O)C(O)(C1CCCC1)C1=CC=CC=C1
 Pathways: Neuroscience
 Receptor: mAChR
 M1
 Boling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

- Haddad E.B., et al. Br J Pharmacol, 1999. 127(2): p. 413-20.

Note: This product is for in vitro research use only