

Nafamostat mesylate

Catalog No: #S2392


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nafamostat mesylate
Brief Description	Inhibitors
Purification	99.76%
Target Name	Serine Protease inhibitor
Calculated MW	539.58
Formulation	C19H17N5O22CH4O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CS(O)(=O)=O.CS(O)(=O)=O.NC(N)=NC1=CC=C(C=C1)C(=O)OC1=CC2=C(C=C1)C=C(C=C2)C(N)=N Pathways: Proteases /Proteasome Receptor: Serine Protease Boiling pt: Melting pt: Solubility: Water: 53 mg/mL (98.22 mM); DMSO: 20 mg/mL (37.06 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sundaram S, et al. Thromb Haemost, 1996, 75(1), 76-82.

Note: This product is for in vitro research use only