

Efavirenz

Catalog No: #S2393


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Efavirenz
Brief Description	Inhibitors
Purification	99.58%
Target Name	HIV Protease inhibitor
Calculated MW	315.68
Formulation	C14H9ClF3NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: C1CC1C#C[C@]2(C3=C(C=CC(=C3)Cl)NC(=O)O2)C(F)(F)F
 Pathways: Microbiology & virology;
 Proteases/Proteasome Receptor: HIV-1 RT
 Boiling pt: Melting pt: Solubility: MEthanol: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Held D.M., et al. J Biol Chem, 2006. 281(35): p. 25712-22.

Note: This product is for in vitro research use only