

## Pixantrone dimaleate

Catalog No: #S2394


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pixantrone dimaleate
Brief Description	Inhibitors
Purification	98.51%
Target Name	Others
Calculated MW	557.21
Formulation	C <sub>25</sub> H <sub>27</sub> N <sub>5</sub> O <sub>10</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer  
 SMILES: C1=CC(=C2C(=O)C3=C(C(=O)C2=C1NCCN)C=CN=C3)NCCN.C(=C\C(=O)O)\C(=O)O.C(=C\C(=O)O)\C(=O)O  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: Blue Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Cavalletti E., et al. Invest New Drugs, 2007. 25(3): p. 187-95.

Note: This product is for in vitro research use only