

Topiroxostat

Catalog No: #S2397


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Topiroxostat
Brief Description	Inhibitors
Purification	98.95%
Target Name	ROS inhibitor
Calculated MW	248.24
Formulation	C13H8N6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:N#CC1=NC=CC(=C1)C1=NNC(=N1)C1=CC=NC=C1Pathways:Immunology/InflammationReceptor:XOBoling pt:Melting

pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sato T, et al. Bioorg Med Chem Lett. 2009 Nov 1;19(21):6225-9.

Note: This product is for in vitro research use only