

## Candesartan Cilexetil

Catalog No: #S2400


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Candesartan Cilexetil
Brief Description	Inhibitors
Purification	99.50%
Target Name	RAAS inhibitor
Calculated MW	610.66
Formulation	C33H34N6O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:CCOC1=NC2=CC=CC(C(=O)OC(C)OC(=O)OC3CCCC3)=C2N1CC1=CC=C(C=C1)C1=CC=CC=C1C1=NNN=N1Pathways:Endocrinology/HormonesReceptor:angiotensin II receptorBoiling pt:Melting pt:Solubility:Ethanol: 4 mg/mL (6.55 mM); DMSO: 122 mg/mL (199.78 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. See S, et al. Am J Health Syst Pharm, 2000, 77(8), 739-746.

Note: This product is for in vitro research use only