

Bromfenac sodium

Catalog No: #S2403


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bromfenac sodium
Brief Description	Inhibitors
Purification	92.00%
Target Name	COX inhibitor
Calculated MW	356.15
Formulation	C15H11BrNNaO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:O.O.O.[Na+].NC1=C(C=CC=C1CC([O-])=O)C(=O)C1=CC=C(Br)C=C1Pathways:Immunology/Inflammation;
NeuroscienceReceptor:COX-2Boiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility ,
please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Shimura M, Yasuda K. Br J Ophthalmol. 2015 Feb;99(2):215-9.

Note: This product is for in vitro research use only