

## Revaprazan hydrochloride


Catalog No: #S2405

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Revaprazan hydrochloride
Brief Description	Inhibitors
Purification	99.91%
Target Name	ATPase inhibitor
Calculated MW	398.9
Formulation	C22H24ClFN4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: CC1N(CCC2=CC=CC=C2)C1=NC(NC2=CC=C(F)C=C2)=NC(C)=C1  
Pathways: Membrane transporter/Ion channel  
Receptor: H<sup>+</sup>/K<sup>+</sup>-ATPase  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Lee JS, et al. J Clin Biochem Nutr. 2012 Sep;51(2):77-83.

Note: This product is for in vitro research use only