

GPR116 Polyclonal Antibody

Catalog No: #40972

Package Size: #40972-1 50ul #40972-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GPR116 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	GPR116 Polyclonal Antibody detects endogenous levels of GPR116 protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GPR116.
Target Name	GPR116
Other Names	GPR116; KIAA0758; Probable G-protein coupled receptor 116
Accession No.	Swiss-Prot: Q8IZF2NCBI Gene ID: 221395
Uniprot	Q8IZF2
GeneID	221395;
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western blot analysis of the lysates from HeLa cells, mouse brain, rat brain using GPR116 antibody.

Note: This product is for in vitro research use only