

Luliconazole

Catalog No: #S2411


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Luliconazole
Brief Description	Inhibitors
Purification	98.27%
Target Name	Antifungal inhibitor
Calculated MW	354.28
Formulation	C ₁₄ H ₉ Cl ₂ N ₃ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1[C@H](S/C(=C(\C#N)/N2C=CN=C2)/S1)C3=C(C=C(C=C3)Cl)Cl Pathways: Microbiology & virology Receptor: Antifungal Boli
 ng pt: Melting pt: Solubility: Ethanol: 5 mg/mL (14.11 mM); DMSO: 71 mg/mL (200.4 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Uchida K, et al. J Infect ChemOthers, 2004, 10(4), 216-219.

Note: This product is for in vitro research use only