


## Description

Product Name	EHop-016
Brief Description	Inhibitors
Purification	99.50%
Target Name	Rho inhibitor
Calculated MW	430.55
Formulation	C25H30N6O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:CCN1C2=CC=CC=C2C2=C1C=CC(NC1=NC(NCCCN3CCOCC3)=NC=C1)=C2Pathways:Cell

Cycle/CheckpointReceptor:Rac1Boling pt:Melting pt:Solubility:DMSO: 86 mg/mL (199.74 mM)Appearance:Off-white Color PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Montalvo-Ortiz BL, et al. J Biol Chem. 2012, 287(16), 13228-13238.

Note: This product is for in vitro research use only