


Description

Product Name	TAS-301
Brief Description	Inhibitors
Purification	98.00%
Target Name	PKC inhibitor
Calculated MW	357.4
Formulation	C23H19NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=CC=C(C=C1)C(=C1C(=O)NC2=CC=CC=C2)C1=CC=C(OC)C=C1Pathways:Cytoskeletal SignalingReceptor:PKCBoling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Muranaka Y, et al. J Pharmacol Exp Ther. 1998 Jun;285(3):1280-1286.

Note: This product is for in vitro research use only