


Description

Product Name	CH 223191
Brief Description	Inhibitors
Purification	98.00%
Target Name	AhR antagonist
Calculated MW	333.39
Formulation	C19H19N5O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: CN1N=CC=C1C(=O)NC1=C(C)C=C(C=C1)N=NC1=CC=CC=C1 Pathways: Metabolism Receptor: AhR Boling pt: Melting pt: Solubility: Ethanol: 4 mg/mL warmed (11.99 mM); DMSO: 66 mg/mL (197.96 mM) Appearance: Orange Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kim SH, et al. Mol Pharmacol. 2006, 69(6), 1871-1878.

Note: This product is for in vitro research use only