

## LAQ824 (Dacinostat)

Catalog No: #S2454


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	LAQ824 (Dacinostat)
Brief Description	Inhibitors
Purification	99.77%
Target Name	HDAC inhibitor
Calculated MW	379.46
Formulation	C22H25N3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer  
 SMILES: C1=CC=C2C(=C1)C(=CN2)CCN(CCO)CC3=CC=C(C=C3)/C=C/C(=O)NO  
 Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair; NF-kB Receptor; HDAC  
 Boiling pt: Melting pt: Solubility: DMSO: 76 mg/mL (200.28 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Atadja P, et al. Y Res, 2004, 64(2), 689-695.

Note: This product is for in vitro research use only