

## IRAK-1/4 Inhibitor

Catalog No: #S2457


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	IRAK-1/4 Inhibitor
Brief Description	Inhibitors
Purification	99.89%
Target Name	Antiviral inhibitor
Calculated MW	395.41
Formulation	C20H21N5O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer SMILES: [O-][N+](=O)C1=CC=CC(=C1)C(=O)NC1=NC2=CC=CC=C2N1CCN1CCOCC1 Pathways: Immunology/Inflammation Receptor: IRAK-1; IRAK-4 Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Powers JP, et al. Bioorg Med Chem Lett. 2006 Jun 1;16(11):2842-5.

Note: This product is for in vitro research use only