


Description

Product Name	RVX-208
Brief Description	Inhibitors
Purification	99.88%
Target Name	Epigenetic Reader Domain inhibitor
Calculated MW	370.4
Formulation	C20H22N2O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:COC1=CC(OC)=C2C(=O)N=C(NC2=C1)C1=CC(C)=C(OCCO)C(C)=C1Pathways:Chromatin/EpigeneticReceptor:BD2Boling
pt:Melting pt:Solubility:DMSO: 74 mg/mL (199.78 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Picaud S, et al. Proc Natl Acad Sci U S A. 2013, 110(49), 19754-19759.

Note: This product is for in vitro research use only