

Rocilinostat (ACY-1215)

Catalog No: #S2489


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rocilinostat (ACY-1215)
Brief Description	Inhibitors
Purification	99.21%
Target Name	HDAC inhibitor
Calculated MW	433.5
Formulation	C24H27N5O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: ONC(=O)CCCCCCNC(=O)C1=CN=C(N=C1)N(C1=CC=CC=C1)C1=CC=CC=C1
 Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair; NF- κ B Receptor: HDAC1; HDAC2; HDAC3; HDAC6; HDAC8
 Boiling pt: Melting pt: Solubility: DMSO: 86 mg/mL (198.38 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Santo L, et al. Blood, 2012, 119(11), 2579-2589.

Note: This product is for in vitro research use only