

MPEP hydrochloride

Catalog No: #S2497


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | MPEP hydrochloride |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | GluR antagonist |
| Calculated MW | 229.7 |
| Formulation | C14H12ClN |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Nervous system
SMILES: Cl.CC1=CC=CC(=N1)C#CC1=CC=CC=C1
Pathways: Neuroscience
Receptor: mGluR
Boiling pt: Melting pt:
Solubility: DMSO: 10 mM
Appearance: Faint Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Inta Dragos, et al. *Neuropharmacology* (2012), 62(5-6), 2034-2039.

Note: This product is for in vitro research use only