

Adjudin (AF-2364)

Catalog No: #S2498


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adjudin (AF-2364)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	335.19
Formulation	C15H12Cl2N4O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: NNC(=O)C1=NN(CC2=C(Cl)C=C(Cl)C=C2)C2=CC=CC=C12 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Mruk DD. Trends Endocrinol Metab. 2008 Mar;19(2):57-64.

Note: This product is for in vitro research use only