

Pritelivir (BAY 57-1293)


Catalog No: #S2504

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pritelivir (BAY 57-1293)
Brief Description	Inhibitors
Purification	98.70%
Target Name	DNA/RNA Synthesis antagonist
Calculated MW	402.49
Formulation	C18H18N4O3S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CN(C(=O)CC1=CC=C(C=C1)C1=CC=CC=N1)C1=NC(C)=C(S1)S(N)(=O)=O
 Pathways: DNA Damage/DNA Repair
 Receptor: helicase primase
 Boiling pt: Melting pt: Solubility: DMSO: 80 mg/mL (198.76 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kleymann G, et al. Nat Med. 2002, 8(4), 392-398.

Note: This product is for in vitro research use only