

Azeliagon

Catalog No: #S2507


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Azeliagon
Brief Description	Inhibitors
Purification	98.00%
Target Name	advanced glycation end products agonist
Calculated MW	532.12
Formulation	C32H38ClN3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CCCCC1=NC(=CN1C1=CC=C(OC2=CC=C(Cl)C=C2)C=C1)C1=CC=C(OCCCN(CC)CC)C=C1 Pathways: Others Receptor: RAGE Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Burstein AH, et al. BMC Neurol. 2014 Jan 15;14:12.

Note: This product is for in vitro research use only