

Galunisertib (LY2157299)

Catalog No: #S2510


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Galunisertib (LY2157299)
Brief Description	Inhibitors
Purification	99.85%
Target Name	TGF-beta/Smad inhibitor
Calculated MW	369.42
Formulation	C22H19N5O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC1=CC=CC(=N1)C1=NN2CCCC2=C1C1=C2C=C(C=CC2=NC=C1)C(N)=O Pathways: Stem

Cells Receptor: TβRI Boling pt: Melting pt: Solubility: DMSO: 74 mg/mL (200.31 mM) Appearance: Faint Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Yingling JM, et al. Proc Am Assoc Y Res. 2006, 47, abstract 250.

Note: This product is for in vitro research use only