

## Mocetinostat (MGCD0103)

Catalog No: #S2512


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Mocetinostat (MGCD0103)
Brief Description	Inhibitors
Purification	99.74%
Target Name	HDAC inhibitor
Calculated MW	396.44
Formulation	C23H20N6O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer SMILES: NC1=CC=CC=C1NC(=O)C1=CC=C(CNC2=NC=CC(=N2)C2=CN=CC=C2)C=C1 Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair; NF-kb Receptor: HDAC1; HDAC11; HDAC2; HDAC3; HDAC4 Boiling pt: Melting pt: Solubility: DMSO: 13 mg/mL (32.79 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Fournel M, et al. Mol Y Ther. 2008, 7(4), 759-768.

Note: This product is for in vitro research use only