

Roxadustat (FG-4592)

Catalog No: #S2515


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Roxadustat (FG-4592)
Brief Description	Inhibitors
Purification	99.50%
Target Name	HIF inhibitor
Calculated MW	352.34
Formulation	C19H16N2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC1=NC(C(=O)NCC(O)=O)=C(O)C2=C1C=C(OC1=CC=CC=C1)C=C2 Pathways: Angiogenesis;

Chromatin/Epigenetic Receptor: HIF-PH Boiling pt: Melting pt: Solubility: DMSO: 70 mg/mL (198.67 mM) Appearance: Faint Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Seeley, TW, et al. U.S. Pat. Appl. Publ. 2007, US 20070292433 A1 20071220.

Note: This product is for in vitro research use only