

GPR52 Polyclonal Antibody

Catalog No: #40983

Package Size: #40983-1 50ul #40983-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR52 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	GPR52 Polyclonal Antibody detects endogenous levels of GPR52 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human GPR52.
Target Name	GPR52
Other Names	GPR52; Probable G-protein coupled receptor 52
Accession No.	Swiss-Prot: Q9Y2T5NCBI Gene ID: 9293
Uniprot	Q9Y2T5
GeneID	9293;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using GPR52 Polyclonal Antibody

Note: This product is for in vitro research use only