

Latanoprost

Catalog No: #S2528


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Latanoprost
Brief Description	Inhibitors
Purification	99.69%
Target Name	GPR agonist; RAR/RXR? antagonist
Calculated MW	432.59
Formulation	C26H40O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC(C)OC(=O)CCC/C=C\C[C@H]1[C@H](C[C@H]([C@@H]1CC[C@H](CCC2=CC=CC=C2)O)O)O

Pathways:Endocrinology/Hormones; MetabolismReceptor:prostanoid selective FP receptor; retinoid X receptor α Boling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Patel SS and CM Spencer Latanoprost. Drugs Aging, 1996. 9(5): p. 363-78.

Note: This product is for in vitro research use only