

Eprosartan Mesylate

Catalog No: #S2531


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eprosartan Mesylate
Brief Description	Inhibitors
Purification	99.93%
Target Name	RAAS inhibitor
Calculated MW	520.62
Formulation	C23H24N2O4SCH4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Endocrine

systemSMILES:CCCC1=NC=C(N1CC2=CC=C(C=C2)C(=O)O)/C=C(\CC3=CC=CS3)/C(=O)O.CS(=O)(=O)O
Pathways: Endocrinology/Hormones
Receptor: AT1 receptor
Boiling pt: Melting pt: Solubility: Ethanol: 4 mg/mL (7.68 mM); DMSO: 104 mg/mL (199.76 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Plosker GL, et al. Drugs, 2000, 60(1), 177-201.

Note: This product is for in vitro research use only