

Tauroursodeoxycholic Acid

Catalog No: #S2532


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tauroursodeoxycholic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	P450
Calculated MW	499.7
Formulation	C26H45NO6S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:C[C@H](CCC(=O)NCCS(=O)(=O)O)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2[C@H](C[C@H]4[C@@]3(CC[C@H](C4)O)C)O)CPathways:MetabolismReceptor:CYP3ABoling pt:Melting pt:173-175CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Paolini M, et al. Hepatology. 1999 Sep;30(3):730-9.

Note: This product is for in vitro research use only