

## Alcaftadine

Catalog No: #S2533


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Alcaftadine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histamine Receptor antagonist
Calculated MW	307.39
Formulation	C19H21N3O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous systemSMILES:CN1CCC(CC1)=C1C2=NC=C(C=O)N2CCC2=CC=CC=C12Pathways:GPCR/G Protein;  
Immunology/Inflammation; NeuroscienceReceptor:H1 receptor; H2 receptor; H4 receptorBoiling pt:Melting pt:Solubility:DMSO:  
40mg/mLAppearance:Faint Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath  
for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Namdar R and C Valdez Alcaftadine. Drugs Today (Barc), 2011. 47(12): p. 883-90.

Note: This product is for in vitro research use only