

## Droxidopa


Catalog No: #S2547

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Droxidopa
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor inhibitor
Calculated MW	213.19
Formulation	C <sub>9</sub> H <sub>11</sub> NO <sub>5</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular System  
 SMILES: C1=CC(=C(C=C1[C@H]([C@@H](C(=O)O)N)O)O)O  
 Pathways: Neuroscience; GPCR/G  
 Protein Receptor: Adrenergic Receptor  
 Boiling pt: Melting pt: 229-232°C  
 Solubility: Ethanol:

## References

1. Mathias CJ. Clin Auton Res, 2008, 18 Suppl 1, 25-29.

Note: This product is for in vitro research use only