

## Succimer

Catalog No: #S2552


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Succimer
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	182.22
Formulation	C <sub>4</sub> H <sub>6</sub> O <sub>4</sub> S <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
SMILES: [C@@H]([C@@H](C(=O)O)S)(C(=O)O)S  
Pathways: Others  
Receptor: Cadmium; Lead; Mercury; Others  
Boiling pt: Melting pt: 196-198°C  
Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Blanusa M, et al. *Curr Med Chem.* 2005;12(23):2771-94.

Note: This product is for in vitro research use only