

## Nandrolone decanoate

Catalog No: #S2556


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nandrolone decanoate
Brief Description	Inhibitors
Purification	99.03%
Target Name	Others
Calculated MW	428.65
Formulation	C28H44O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:CCCCCCCCC(=O)O[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CCC4=CC(=O)CC[C@H]34)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:33-37CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Pozzi R, et al. Toxicol Mech Methods. 2013 May;23(4):289-96.

Note: This product is for in vitro research use only