

## Allylestrenol

Catalog No: #S2559


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Allylestrenol
Brief Description	Inhibitors
Purification	99.50%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	300.48
Formulation	C21H32O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine system  
SMILES:CC12CCC3C(CCC4=CCCC34)C1CCC2(O)CC=C  
Pathways:Endocrinology/Hormones  
Receptor:ER  
Boiling pt:Melting pt:  
Solubility:DMSO: 10 mM  
Appearance:White Crystal  
Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
Stock solution can be stored below -20 °C for several month.

## References

1. Madauss KP, et al. Med Res Rev. 2007 May;27(3):374-400.

Note: This product is for in vitro research use only