

O-Acetyl-L-carnitine hydrochloride

Catalog No: #S2563


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	O-Acetyl-L-carnitine hydrochloride
Brief Description	Inhibitors
Purification	96.99%
Target Name	Others
Calculated MW	239.7
Formulation	C ₉ H ₁₈ ClNO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: CC(=O)O[C@H](CC(=O)O)C[N+](C)(C)C.[Cl-]
 Pathways: Others
 Receptor: CRAT
 Boiling pt: Melting pt: 194°C
 Solubility: DMSO: 10 mM
 Appearance: White Paste solid
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only