

Etonogestrel

Catalog No: #S2567


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Etonogestrel
Brief Description	Inhibitors
Purification	99.59%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	324.46
Formulation	C22H28O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC[C@]12CC(=C)[C@H]3[C@H]([C@@H]1CC[C@]2(C#C)O)CCC4=CC(=O)CC[C@H]34Pathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:182-184CSolubility:DMSO: 10 mMAppearance:Transparent CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. McNicholas C, et al. Obstet Gynecol. 2015 Mar;125(3):599-604.

Note: This product is for in vitro research use only