

Felypressin

Catalog No: #S2568


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Felypressin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Vasopressin Receptor agonist
Calculated MW	1040.22
Formulation	C46H65N13O11S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C1C[C@H](N(C1)C(=O)[C@@H]2C[SSC][C@@H](C(=O)N[C@H](C(=O)N[C@H](C(=O)NC(C(=O)N[C@H](C(=O)N2)CC(=O)N)CCC(=O)N)CC3=CC=CC=C3)CC4=CC=CC=C4)N)C(=O)N[C@H](CCCCN)C(=O)NCC(=O)NPathways:GPCR/G ProteinReceptor:vasopressin receptor
 1aBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White Cotton PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Akkerman JW, van den Hoven WE. Drugs. 1987;33 Suppl 2:53-9.

Note: This product is for in vitro research use only