

Chlorotrianisene

Catalog No: #S2569


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chlorotrianisene
Brief Description	Inhibitors
Purification	99.89%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	380.86
Formulation	C23H21ClO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:COC1=CC=C(C=C1)C(Cl)=C(C1=CC=C(OC)C=C1)C1=CC=C(OC)C=C1Pathways:Endocrinology/HormonesReceptor:ERBoiling

pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kupfer D, Bulger WH. FEBS Lett. 1990 Feb 12;261(1):59-62.

Note: This product is for in vitro research use only