

Benactyzine hydrochloride

Catalog No: #S2571


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Benactyzine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR antagonist
Calculated MW	363.88
Formulation	C20H26ClNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: Cl.CCN(CC)CCOC(=O)C(O)(C1=CC=CC=C1)C1=CC=CC=C1 Pathways: Neuroscience Receptor: mAChR Boling pt: Melting

pt: 177-179°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Amitai G, et al. Mol Pharmacol. 1987 Nov;32(5):678-85.

Note: This product is for in vitro research use only