

Guanethidine Monosulfate

Catalog No: #S2575


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Guanethidine Monosulfate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Norepinephrine inhibitor
Calculated MW	296.39
Formulation	C10H24N4O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine system SMILES:OS(O)(=O)=O.NC(N)=NCCN1CCCCC1 Pathways:Neuroscience Receptor:Norepinephrine (NE) Boiling pt:Melting pt:Solubility:Soluble in Water Appearance:White Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fabiani ME, et al. Pharmacol Res. 1996 Mar;33(3):171-80.

Note: This product is for in vitro research use only