

Dienogest


Catalog No: #S2577

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dienogest
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	311.42
Formulation	C ₂₀ H ₂₅ NO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CCC3=C4CCC(=O)C=C4CC[C@H]3[C@@H]1CC[C@]2(CC#N)O

Pathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:210-214CSolubility:DMSO: 62 mg/mL (199.08 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sasagawa S, et al. Steroids, 2008, 73(2), 222-231.

Note: This product is for in vitro research use only