

## Reboxetine

Catalog No: #S2579


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Reboxetine
Brief Description	Inhibitors
Purification	99.72%
Target Name	Norepinephrine inhibitor
Calculated MW	313.39
Formulation	C19H23NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous

system SMILES: CCOC1=CC=CC=C1OC(C2=CC=CC=C2)C3=CC=CC=C3 Pathways: Neuroscience Receptor: Norepinephrine (NE) Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Rauhut AS, et al. J Pharmacol Exp Ther. 2002 Nov;303(2):664-72.

Note: This product is for in vitro research use only